

Clinical Photos



Before



After



Before



After



Before




After




Before



After



Before



After

For More Clinical Photos



1. Scan the QR code.
2. Select SHINY PLUS in 'Search By Device' menu.

The Safest Treatment Micro Pulsed Light

Indications

Pigmented Lesions

Solar lentigines | Freckles | Melasma
Hyperpigmentation (PIH)

Rejuvenation

Wrinkles | Enlarged pores

Acne

Acne Vulgaris

Hair Removal

Leg hair | Axillary hair | Bikini hair

Vascular Lesions

Telangiectasia | Flushing | Rosacea

Want to Know More?

Want to know more? Just scan the QR code!



www.ids.co.kr



Tx Video

ids



SHINY PLUS

The Safest Treatment
by Micro Pulsed Technology


FDA CE

SHINY PLUS


SHINY PLUS is an advanced MPL (Micro Pulsed Light) device designed to treat a variety of skin concerns with high efficacy and minimal side effects. Its uniquely engineered sapphire tip, featuring 140 micro-cells, helps prevent common IPL issues like “tiger marks.”

Micro Pulsed Technology

The multi-cell sapphire design isolates light energy within micro-cells, leaving safety zones between each. This design protects healthy skin, making it safe for all skin types.



Excessive IPL can cause large blisters, risking PIH or scarring. SHINY PLUS minimizes this by producing tiny vesicles instead, reducing the chance of serious side effects.



Immediate

After 24 hours

48 hours after an intentional excessive SHINY PLUS treatment, only small vesicles appeared—instead of a single large blister.

Seven Filters

Seven interchangeable filters allow precise, customized treatments for various skin concerns and all skin types



Frequently Asked Questions

Q: Is SHINY PLUS safe for all skin types?

A: Absolutely. With its advanced Micro Pulsed Technology, multi-cell sapphire tip, and seven specialized filters, SHINY PLUS ensures the safest treatments for all skin tones—from light to dark.

Q: What are the benefits of the 7 filters?

A: The seven filters allow tailored treatments for various skin conditions and optimize energy delivery for safe, targeted results.

Q: Does the treatment hurt?

A: Most patients find the treatment comfortable. The sapphire tip provides a cooling effect that enhances comfort during sessions.

Q: How long does a treatment session take, and how many sessions are needed?

A: Treatment time varies depending on the area and concern but typically lasts between 15 to 30 minutes. While visible results can appear after just one session, a series of 3–6 treatments is usually recommended for optimal and lasting results.

Q: Is there any downtime after SHINY PLUS treatment?

A: There is little to no downtime. Some redness may occur, but it usually subsides within a few hours.

