

CRYO-ZET

Air Cooling System

CRYO-ZET is an air cooling device designed to be used with lasers during dermatological procedures. It minimizes skin damage through compatibility with various lasers, reduces pain and recovery time, and ensures even skin protection through its advanced cooling capabilities.

With the capability to chill the air to a maximum of -35 degrees Celsius, it delivers a robust cooling effect both during and after laser treatments.

Notably, CRYO-ZET operates without the need for consumables, ensuring cost-effectiveness, simplicity, and convenience in its usage.



For international use only



Feature

Protect skin

Minimizes skin damage and reduces pain and recovery time.

Air cooling

Utilizes an airflow control system to regulate temperature.

Temperature regulation

Achieves cooling air temperatures down to a maximum of -35 degrees Celsius.

Quality compressor

Employs a high-quality compressor for rapid cooling.

Automation

Incorporates an auto-adjustable control system and auto-detecting system to maintain optimal conditions.



Technical specification

Air Level	10 step
Air Flow Output	1000L/min
Temperature	0 to -35 °C
AC Input	220-240 V, 50/60 Hz
Dimension (W*L*H)	400 mm (W) X 610 mm (L) X 963 mm (H)
Weight	50 Kg