

Indication

Diamond Peel

Mechanical Exfoliation

- Removes dead skin cells and impurities
- Refines skin texture
- Minimizes pores
- Promotes healthy cell turnover
- Enhances skin brightness
- Prepares the skin for better absorption of active ingredients

Ultrasonic

Sonophoresis & Micro-vibration Massage

- Enhances the delivery of active ingredients into the skin
- Stimulates blood circulation and promotes lymphatic drainage
- Improves skin tone, elasticity, and absorption of skincare products

Hydrofacial

Hydration & Suction-based Cleansing

- Purifies the skin
- Replenishes moisture
- Removes dead skin cells, excess sebum, and impurities
- Improves skin clarity, texture, and overall radiance

Radiofrequency

Collagen Production & Mild Lifting

- Stimulates collagen and elastin production
- Promotes skin tightening and improves overall elasticity
- Reduces fine lines, mild sagging, and early signs of aging

Liner HiFU*

Lifting & Firming

- Achieves a natural lifting effect
- Tightens and firms sagging skin

* Optional

Technical specification

| | | |
|-------------------|---------------------------------------|---|
| Diamond Peel HP | #140 / #170 | |
| Hydrofacial HP | Regular / Narrow | |
| Ultrasonic HP | 1 MHz / 3 MHz | |
| Radiofrequency HP | 1 MHz | |
| Liner HiFU HP* | Standard Line: 3.0, 4.5 mm | Optional Line: 2.0, 6.0, 8.0 mm Dot: 1.5, 3.0, 4.5, 6.0, 9.0, 10.0, 13.0 mm |
| Vacuum Pump | Up to 650 mmHg (±10 %) | |
| Power Consumption | 80 W | |
| Aqua Waste Tank | 1200 cc | |
| Input Voltage | 100/220 VAC, 50/60 Hz | |
| Dimension (W*L*H) | 460 mm (W) X 420 mm (L) X 1100 mm (H) | |
| Weight | 30 Kg | |

* Optional

HYDRO PLUS

All-in-One Skin Treatment Platform



Diamond Peel Handpiece



Hydrofacial Handpiece



Ultrasonic Handpiece



Radiofrequency Handpiece



Liner HiFU Handpiece*

* Optional

HYDRO PLUS, a next-generation all-in-one skin treatment platform, is engineered to elevate facial rejuvenation through a synergistic combination of technologies. Designed to target multiple layers of the skin, HYDRO PLUS integrates five specialized handpieces—Diamond Peel HP, Hydrofacial HP, Ultrasonic HP, Radiofrequency HP, and Liner HiFU HP—for comprehensive care from exfoliation to lifting.

This versatile system addresses a wide range of skin concerns including dullness, dehydration, loss of elasticity, and signs of aging. By combining mechanical peeling, deep serum infusion, and energy-based tightening, HYDRO PLUS empowers practitioners to deliver visible, long-lasting results tailored to individual skin needs across all skin types.



Diamond Peel for Mechanical Exfoliation

Diamond Peel offers precise mechanical exfoliation using a diamond-coated tip to gently remove dead skin cells and impurities from the skin surface. This non-invasive procedure helps refine skin texture, minimize pores, and promote healthy cell turnover. It enhances skin brightness and prepares the skin for better absorption of active ingredients in subsequent treatments. Ideal for all skin types, Diamond Peel leaves the skin smoother, clearer, and more radiant.

Available Configurations

- #140
- #170



Diamond Peel Handpiece

Hydrofacial for Hydration and Suction-based Cleansing

Hydrofacial combines hydration therapy with suction-based cleansing to deeply purify the skin while replenishing moisture. Its unique spiral tip design effectively removes dead skin cells, excess sebum, and impurities while delivering moisture deep into the skin. This dual-action process improves skin clarity, texture, and overall radiance. Suitable for all skin types, Hydrofacial leaves the complexion feeling clean, smooth, and deeply hydrated.

Available Configurations

- Regular (for face)
- Narrow (for nose)



Hydrofacial Handpiece

Ultrasonic for Sonophoresis and Micro-vibration Massage

Ultrasonic uses high-frequency sound waves to enhance the delivery of active ingredients into the skin through a process called sonophoresis. Simultaneously, gentle micro-vibration massage stimulates blood circulation and promotes lymphatic drainage. This dual-action treatment helps improve skin tone, elasticity, and absorption of skincare products. It is ideal for revitalizing dull, tired skin and achieving a more radiant complexion.

Available Configurations

- 1 MHz
- 3 MHz



Ultrasonic Handpiece

| | | |
|--|--|--|
| Diamond Peel HP Mechanical exfoliation | Hydrofacial HP Hydration, Suction-based cleansing | Ultrasonic HP Sonophoresis, Micro-vibration massage |
| Radiofrequency HP Collagen production, Mild lifting | Liner HiFU HP* Lifting, Firming | |

* Optional

Radiofrequency for Collagen Production and Mild Lifting

Radiofrequency (RF) delivers controlled thermal energy deep into the skin to stimulate collagen and elastin production. This non-invasive treatment promotes skin tightening and improves overall elasticity, resulting in a firmer and more lifted appearance. It is especially effective for reducing fine lines, mild sagging, and early signs of aging. Safe for most skin types, RF provides gradual, natural-looking rejuvenation with no downtime.

Available Configurations

- 1 MHz



Radiofrequency Handpiece

Liner HiFU for Lifting and Firming

HiFU targets the deeper layers of the skin, including the SMAS layer, to stimulate collagen remodeling and achieve a natural lifting effect. As a result, it tightens and firms sagging skin without the need for invasive procedures.

Especially, HYDRO PLUS is equipped with Liner HiFU cartridges that deliver line-type HiFU energy, providing a more comfortable and pain-free experience compared to traditional dot-type HiFU systems.

Available Configurations

- | | |
|---------------------------------|---|
| Standard • Line: 3.0, 4.5 mm | Optional • Line: 2.0, 6.0, 8.0 mm • Dot: 1.5, 3.0, 4.5, 6.0, 9.0, 10.0, 13.0 mm |
|---------------------------------|---|



Liner HiFU Handpiece*

* Optional